

Use of the Axiocam 208 on the Leica S8APO in the insect lab

(A) For illustration purposes on wall-mounted screen and simple photo & video applications

Note 1: Images/videos are saved to the memory stick on USB hub. Please clean the stick up after use.

Note 2: The keyboard attached to the USB hub is Swiss, but the key assignment is American.

- 1) Turn on wall-mounted monitor (indented button at the back, bottom right).
- 2) Turn on camera with switch on the black power chord.
- 3) Turn on light on the microscope stand (try different lighting options for best results).
- 4) Set microscope oculars to zero correction for matching focus of microscope and camera.
- 5) Call up the on screen menu (cogwheel button on camera), make sure Exposure and White Balance are on 'Auto'.
- 6) Put your object under the microscope and adjust focus and/or magnification.
Now you should have an image on screen. For simple observation that should be fine.
For photography or filming, you may want to make adjustments (see below).
- 7) To take a picture, click 'Snap', to record a movie, click 'Record'.

Basic adjustments under 'Home' menu:

- 8) Brightness: Set Exposure to 'Manual' & adjust exposure time with slider or by typing a value.
- 9) Color: Set White Balance to 'Manual' & adjust with slider (left = warmer colors, right = colder colors). 'Push' button will reset to default.

Advanced adjustments under 'Settings' -> Camera

- 10) Gamma (a correction function to better match computer-generated colors to colors perceived by the human eye). The default setting (0.45) works well, but can be adjusted according to taste.
- 11) Optimize options: I recommend to have 'Sharpen' ON, the rest OFF. But check it out.
- 12) Image Orientation: Is self-explanatory.
- 13) Scale Bar: Is only useful when Camera is operated via PC & ZEN core software.

Advanced adjustments under 'Settings' -> Microscope

- 14) Shading Correction: If necessary, follow on-screen instructions. Probably not needed in most cases.

Advanced adjustments under 'Settings' -> Operating System

- 15) Language: Please leave on English.
- 16) Date & Time: Please leave (unless summer time/winter time adjustment is necessary).
- 17) File Options: Custom fields allow you to set a prefix for all pictures you take in a session (e.g. your name). The automatic addition of date & time is useful, but can be turned off.
For 'File Type' use JPEG for standard imaging and smaller files, TIFF if you want to do a lot of further image editing (generates larger files).
- 18) HDMI resolution: Leave on 4K (monitor can handle it). May need to be changed if camera were to be attached to other monitor (not foreseen).
- 19) Bright Pixel Correction: Do not use.
- 20) Firmware Update: Ignore.

(B) For operating the camera from your laptop computer

Install ZEN core software on your laptop (current version is 3.5, stored on memory stick in the USB hub & downloadable: www.mikroskoptechnik.ch/download/zeiss/ZEISS_ZEN_Core_V3.5-MTD.zip).

I also put it into our folder on our server:

Q:\Abteilungsprojekte\eco\ResGroups\Vorburger\Zeiss_Axiocam_Software\ZEISS_ZEN_Core_V3.5

Instructions for offline installation:

<https://mikroskoptechnik.ch/de/kundenbereich/zeisssoftware/coreoffline/index.php>

Connect camera to your laptop using the labelled cable provided at the microscope table.

Self-educate on using the camera from your computer. A good starting point is the “Kurzanleitung” (in German) provided by Diethelm Microscopes, available on our server space:

Q:\Abteilungsprojekte\eco\ResGroups\Vorburger\Zeiss_Axiocam_Software\ZEISS_ZEN_Core_V3.5

An English Quick Guide seems to be available here:

file:///C:/Users/vorburch/Downloads/EN_quick-guide_ZEN-blue-edition_first-steps.pdf

There may be other online resources I have not discovered yet.

Note: We use the free version of the ZEN software. Therefore, not all functions of the software are available.

CV, 18/1/2023