

# Pea Aphid Microsatellite Multiplex

(7 Loci: AIA09M, AIB07M, AIB08M, AIB12M, ApF08M, ApH08M, ApH10M from Caillaud et al. 2004)

## Material:

- PCR plates/strips
- ABI plates
- 1.5ml Tubes
- Réservoir
- Cold elements/ cool boy
- extra waste bag for HiDi waste

## Chemical:

- Qiagen Master-Mix
- Primers
- HiDi: **Suspected of damaging the unborn child, Suspected of causing cancer**  
**May cause damage to organs through prolonged or repeated exposure**  
Work with lab coat and with nitrile gloves and the best would be to work under the hood.

## PCR:

- Fill out *MP\_Pea\_aphid* excel sheet and pipette mastermix accordingly.
- Distribute 10 ul mastermix to each well of the PCR plate/strip, then add DNA (1 ul per well)
- Centrifuge everything down
- Let PCR run with PCR-program: *MP\_pea\_aphid*. If this is not saved on the machine yet, program the machine according to the *MP\_Pea\_aphid* excel sheet and save program in the Vorbürger lab folder.

## ABI:

- tell Esther or Katri in advance that you are going to prepare a ABI plate!!
- dilute PCR product 1:20 -> 1ul PCR + 19ul autoclaved MiliQ
- centrifuge plate quickly
- **Work now under the hood or put everything on a cold element and wear lab coat + nitrile gloves!!**
- put the ABI plate in a cool boy to cool down
- Prepare HiDi+Liz500 mix: - 9.9 ul HiDi \*(Sample number + 10% sample number) + 0.1 ul Liz500\*(Sample number + 10% sample number) **Hidi waste goes in separate waste bag!**
- Pipette in each well of an ABI plate 10ul of HiDi+Liz500 mix
- add 0.5 ul of diluted PCR product
- If one of the even numbered rows is empty, pipette 10ul of HiDi in each well of this row
- Fill out ABI data sheet and save it as a text and excel file under [Q:\Querprojekte\MolGen\ABI new 2010\data sheets](#)
- Put the plate in the fridge of the ABI room
- Put the waste bag under the hood and throw it away after one week
- Tell Esther or Katri that the plate is ready to let run