

## Spiroplasma MLST

### 1 DIAGNOSTIC PCR

Pre-heat lid to 95°C

3 min at 95°C

10x cycle of

|         |          |
|---------|----------|
| 95°C    | for 30 s |
| 65-56°C | for 30 s |
| 72°C    | for 60 s |

25x cycle of

|      |          |
|------|----------|
| 95°C | for 30 s |
| 55°C | for 30 s |
| 72°C | for 60 s |

6 min at 72°C

Hold at 10°C

Hot start (waiting for the block to reach temperature) is not required, ramp rate 3-5°C/sec

➔ For RpoBF1.ixod / RpoBR1: set elongation time to 120 s in both cycles

### 2 REAGENTS PER REACTION

5 µl Promega M7822 GoTaq G2 Green Master Mix

2 µl Nuclease free water

1 µl Forward primer at 10 µM

1 µl Reverse primer at 10 µM

1 µl template DNA

### 3 PRIMERS

#### 3.1 16S

Together, the 16S primers cover regions 1-8 out of 9 variable regions in Spiro 16S

| Name      | OligoID & length | Sequence                          | Comments               |
|-----------|------------------|-----------------------------------|------------------------|
| 23F       | 2366024<br>23 bp | 5' CTC AGG ATK AAC GCT GGC GGC AT |                        |
| RiboTKSSR | 2366025<br>20 bp | 5' TAG CCG TGG CTT TCT GGT AA     | Product size: ~ 481 bp |
| 16STF1    | 2366026<br>23 bp | 5' GGT CTT CGG ATT GTA AAR GTC TG |                        |
| 16STR1    | 2366027<br>20 bp | 5' GGT GTG TAC AAG ACC GAG AA     | Product size: ~ 975 bp |

### 3.2 ITS

These are improvements of primers SP-ITS-JO4 (5' GCC AGA AGT CGG TGT CCT AAC CG) and SP-ITS-N55 (5' ATT CCA AGG CAT CCA CCA TAC M). ITS was not used to calculate phylogenies, as the primer amplified products of different lengths in some samples and the sequencing data was not useable.

| Name       | OligoID & length | Sequence                         | Comments             |
|------------|------------------|----------------------------------|----------------------|
| SP-ITS-CVf | 2371743<br>22 bp | 5' GAT GGA TCA CCT CCT TTC TAT G |                      |
| SP-ITS-CVr | 2371744<br>19 bp | 5' CAG CTC CCC GAA SCT TAT C     | Product size: 296 bp |

### 3.3 RpoB

RpoBF1.ixod and RpoR1 produce low yield (weak bands on TBE gel) for some Spiroplasma-positive pea aphid samples and also caused problems during sequencing in the past. Always run RpoBF1.ixod and RpoBR1 with an elongation time of 120s. Product of the two primers overlap.

| Name         | OligoID & length | Sequence                                 | Comments  |
|--------------|------------------|--|---|
| RpoB.cv.F2   | 2375067<br>20 bp | 5' TCC ATC ACG TGG AAC TTG AC            |   |
| RpoB.cv.R2   | 2375068<br>23 bp | 5' CGG TAA GGT GTT TCA ATA AAA CC        | Product size 1'179 bp                                   |
| RpoBF1.ixod* | 2375069<br>26 bp | 5' ATG GAT CAA ACT AAT CCA TTA<br>GCA GA | Requires 120 s elongation step                          |
| RpoBR1*      | 2375070<br>23 bp | 5' CCA AAT CTT TGT CCA CCA TTT TG        | Requires 120 s elongation step<br>Product size 1'784 bp |

### 3.4 ApDNAa

| Name      | OligoID & length | Sequence                             | Comments            |
|-----------|------------------|--------------------------------------|---------------------|
| ApDnaA_F1 | 2429956<br>23 bp | 5' ATT CTT CAG TAA AAA TGC TTG<br>GA |                     |
| ApDnaA_R1 | 2429957<br>23 bp | 5' ACA CAT TTA CTT CAT GCT ATT<br>GA | Product size 447 bp |